

The Lombart **CS-IVs** Chair & Stand...

Made in America, the **CS-IVs** was developed exclusively for Lombart. The **CS-IVs Chair and Stand** features streamlined styling for today's modern office, at an economical price point.

The Chair

Richly upholstered and ergonomically designed, the **Signature Series** manual reclining chair with pneumatic assist lays out flat for a full 180 degree table effect. Retractable armrests and footrest facilitate access while the articulating headrest is easily positioned and securely locked. The chair rotates a full 360 degrees and locks in any position. In addition, the optional footswitch can be conveniently positioned for ease of elevation.


The Stand

Streamlined yet stable, the **CS-IVs Stand** provides a solid foundation for your instruments. The **counter-balanced** refractor and slit lamp arms allow for quick and effortless positioning and can be secured with a single locking lever. Optionally available are your choice of fixed or counter-balanced keratometer arms. The stylish and flexible overhead lamp can be optionally positioned and its intensity is easily adjusted via the console mounted rheostat.

The Modular Console

The heart of the **CS-IVs Stand** is the **ETL (Electronic Testing Lab) Listed modular console** with features you would expect to find on more expensive units. Included are three recharging wells for hand instruments, chair elevation and auto-return switches and two rheostats for adjustment of the binocular indirect ophthalmoscope and overhead lamp

intensity. We've even included a chart projector on/off switch.




OPTIONAL:

Fixed K-Arm or Counter-Balanced K-Arm



CS-IVs with optional Workstation



CS-IVs with optional fixed K-arm