

710

Minor procedures  
chair/table with  
pillow headrest



**Rehance®**

# 710



Features foldaway headrest and reclines 180° for minor procedures

## Surgical Headrest

A non-articulating pillow headrest is standard. An optional special articulating pillow headrest features a sliding, removable magnetic pad and a fold-away design for otology. Available on Model 710 only.

## Powered Base And Articulating Top

Use finger switches and corded foot switch to raise, lower or recline chair up to 180° without breaking contact with patient. Eliminates wasted movements. An auto return button automatically returns patient to the low, upright position.

## Memory Positioning/Safety Switch

Standard equipment includes both a memory switch and a safety switch. The memory switch enables you to pre-select the chair's most used position, while the patented programmable safety switch prevents children from operating the chair by disabling side switches and foot switches.

## Membrane Switches

Membrane switches enhance the Model 710's sleek design and contribute to its durability and safety. They last 10-times longer than their mechanical switches, eliminate the need for high voltage, and allow for limited programming of chair functions. Back lighting eases work in a dimmed room. Mechanical switches available as option.

## Dual Rotation Locks

Easily accessible levers located on both sides of seat enable practitioner to swivel chair 330° without changing position.

## Large, Hinged Footstep

The 8" x 13" footstep, which is padded underneath for comfort, supports patient's legs when seated, and provides safe ingress and egress. Tilts out of way. An upholstered pad supports the feet of patients for standard procedures.

## Armrests

Armrests lock in place. Release mechanism allows convenient egress and ingress from either side of chair. Padded underneath for comfort.

## Upholstery

Match office or treatment room decor with our selection of 14 easy-to-maintain vinyl colors. Other colors can be supplied at additional cost. Available with either smooth or the designer-style cushions shown on cover. Base is available in two standard enamel colors: #17 light beige and #15 silver smoke. All materials are fire retardant. To request a color selection card, contact your dealer or call Reliance at 1-800-735-0357.

## The UL Mark Of Quality

The Model 710 chair has passed rigorous testing by Underwriters Laboratories, which assures its quality and safety for you and your patients.

## Dimensions and Specifications

Base dimensions 22 1/8" x 24" diameter. Overall width: 24 1/4". Overall depth in chair position (excluding headrest): 40". Table position: 64". Seat: 18 1/2" (wide) x 18" (deep). Distance between armrests: 19". Choice of two base heights: HIGH: 22 1/2" lowered, 34" raised. LOW: 20" lowered, 29" raised. Shipping weight: ca. 290 lbs. Electrical: 115 volts, 60 Hz, 8.0 amperes.

**HS** RELIANCE  
INTERNATIONAL

Precision by Tradition.